



# ELK Biotechnology

RFP-Tag Rabbit pAb

Catalog NO.: EA008

For research use only.

## Overview

---

Product name	RFP-Tag Rabbit polyclonal antibody
Source	Rabbit
Applications	<b>WB</b>
Species reactivity	<b>N/A</b>
Recommended dilutions	<b>WesternBlot:1/5000</b> <b>NOTE: Optimal dilutions should be determined by the end user.</b>
Immunogen	Recombinant Protein
Species	N/A
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>N/A</b>
GeneID (Human)	N/A
Human Swiss-Prot No.	N/A
Cellular localization	N/A
Alternative Names	N/A
Background	The red fluorescent protein cloned from Discosoma coral (DsRed or drFP583) holds great promise for biotechnology and cell biology as a spectrally distinct companion or substitute for the green fluorescent protein(GFP) from the Aequorea jellyfish.

